


Closed circuit drainage bags



Description

- PVC sterile drain bag, 5000 ml capacity
- Non return valve
- “Pasteur” type drip chamber with air filter
- Wide diameter connection tube with anti-kneeling spiral
- End conic connector with protective cap and needle free sample port
- “T” cross outlet
- Bag hanger
- Disposable, sterile, double packaged

Advantages

- Closed circuit assures safe drainage reducing the risk of urinary tract infections
- “Pasteur” type drip chamber prevents intraluminal migration of infections
- Needle-free sample port
- High quality “T” valve assures leakage-free bag
- Non-return valve for safe use
- Sterile double packaging for aseptic use

Code	Description	Capacity	Sale box
MS/00060G	Closed circuit sterile drain bag	5000 ml.	100

