


CHIMED SISTEMI MEDICALI
ITALIA

UROCHIM PRODUCT LINE

Urine collection bags



Index

Company profile				3
Product description – Advantages				4
Euro bag	2000ml	without outlet	5
Eurobasic bag	2000ml	without outlet	5
Europlus bag	2000ml	push-pull valve	6
Europlus T-tap valve bag	2000ml	T-tap valve	7
Europlus + bag	2000ml	T-tap valve	8
“Personal” leg bag	750ml	push-pull and T-tap valve	9
Euro 5 bag	5000ml	T-tap valve	10
Euro 5 emptying bag	5000ml	without outlet	11

COMPANY PROFILE

With its 40-year experience in the biomedical field, Chimed S.r.l. manufactures medical devices adhering to the highest quality standards.



Highest consideration is given to the components used in the production cycle, with the aim of ensuring the final products meet the expectations of customers and are compliant with the applicable regulations.



The production process is certified according to ISO 13485:2016 regulations and takes place in Chimed's own clean room, employing modern equipment specifically designed for the medical sector. The procedures in use guarantee

the highest quality of finished products in order to provide the final user with the utmost degree of safety and reliability.

The collection bags manufactured by Chimed have been projected to offer a high practicality and an easy handling, ensuring at the same time a safe and effective use. The various models available meet the multiple needs of patients in both homecare and hospital settings.




Product description

- PVC urine collection bags 750, 2000, or 5000ml capacity, sterile and non-sterile.
- Available with or without non-return valve
- Non-adherent inner surface
- Space on front of bag to indicate patient's details
- Pre-punched hanger holes for attachment to bed support
- PVC connecting tube in various lengths, with catheter-cone connector and protective cap
- With or without drain valve
- Push-pull or T-tap drain valve
- Sterile version packed in blister


Strong points

- High-quality drain valves for a leakage-risk-free product
- Anti-kinking PVC connecting tube
- Devices manufactured using latex-free and phthalates-free materials
- Universal catheter-cone connector for the insertion of anti-slip bladder catheters
- High-frequency double perimeter welding
- Bags entirely manufactured and packaged in clean room
- The non-return valve avoids urine backflow and helps reducing the risk of infections
- The double packaging allows a correct handling of the sterile product
- Devices manufactured at Chimed's premises in Livorno, Italy in compliance with all the applicable regulations (ISO 13485 certified company)


Euro bag 2000ml without outlet

REF	Description	Capacity	Sales box	Picture
101100	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 90cm, non-return valve. Non sterile.	2000ml	250 bags	
101104	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 90cm, non-return valve. Non sterile. Pack of 30 bags.	2000ml	6 packs of 30 bags	
101140	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 130cm, non-return valve. Non sterile.	2000ml	200 bags	
101144	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 130cm, non-return valve. Non sterile. Pack of 30 bags.	2000ml	5 packs of 30 bags	


Eurobasic bag 2000ml without outlet

REF	Description	Capacity	Sales box	Picture
101200	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 90cm. Non sterile.	2000ml	250 bags	
101204	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 90cm. Non sterile. Pack of 30 bags.	2000ml	8 packs of 30 bags	
101240	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 130cm. Non sterile.	2000ml	200 bags	
101244	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 130cm. Non sterile. Pack of 30 bags.	2000ml	7 packs of 30 bags	




Europlus bag 2000ml push-pull valve

REF	Description	Capacity	Sales box	Picture
102100	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 90cm, with push-pull valve, non-return valve. Non sterile.	2000ml	200 bags	
102104	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 90cm, with push-pull valve, non-return valve. Non sterile. Pack of 30 bags.	2000ml	6 packs of 30 bags	
102140	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 130cm, with push-pull valve, non-return valve. Non sterile.	2000ml	180 bags	
102144	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 130cm, with push-pull valve, non-return valve. Non sterile. Pack of 30 bags.	2000ml	6 packs of 30 bags	
112100	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 90cm, with push-pull valve, non-return valve. Sterile.	2000ml	100 bags individually packaged	
112144	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 130cm, with push-pull valve, non-return valve. Sterile.	2000ml	100 bags individually packaged	





Europlus T-tap valve bag 2000ml T-tap valve

REF	Description	Capacity	Sales box	Picture
102200	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 90cm, with T-tap drain valve, non-return valve. Non sterile.	2000ml	200 bags	
102204	Urine collection bag 2000ml capacity with connecting tube, Ø 4.8x6.8mm, length 90cm, with T-tap drain valve, non-return valve. Non sterile. Pack of 30 bags.	2000ml	6 packs of 30 bags	
102240	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 130cm, with T-tap drain valve, non-return valve. Non sterile.	2000ml	180 bags	
102244	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 130cm, with T-tap drain valve, non-return valve. Non sterile. Pack of 30 bags.	2000ml	5 packs of 30 bags	
112200	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 90cm, with T-tap drain valve, non-return valve. Sterile.	2000ml	100 bags individually packaged	
112240	Urine collection bag 2000ml capacity with connecting tube Ø 4.8x6.8mm, length 130cm, with T-tap drain valve, non-return valve. Sterile.	2000ml	100 bags individually packaged	

Europlus+ bag 2000ml T-tap valve



REF	Description	Capacity	Sales box	Picture
104143	Urine collection bag 2000ml capacity with connecting tube Ø 6x9mm, length 130cm, open jaw blue slide clamp, needle-free sampling port and T-tap drain valve. Non sterile.	2000ml	100 bags	
114143	Urine collection bag 2000ml capacity with connecting tube Ø 6x9mm, length 130cm, open jaw blue slide clamp, needle-free sampling port and T-tap drain valve. Sterile.	2000ml	100 bags individually packaged	
114147	Urine collection bag 2000ml capacity with connecting tube Ø 6x9mm, length 130cm, open jaw blue slide clamp, needle-free sampling port and T-tap drain valve. Equipped with bag hanger. Sterile	2000ml	100 bags	

Personal leg bag 750ml push-pull valve and T-tap valve



REF	Description	Capacity	Sales box	Picture
103104	Urine collection leg bag 750ml capacity with tube Ø 4.8x6.8mm, length 10cm, push-pull valve. Non sterile. Pack of 10 bags. <i>(one pack of 10 pieces contains two straps with button)</i>	750ml	40 packs of 10 bags	
103144	Urine collection leg bag 750ml capacity with tube Ø 4.8x6.8mm, length 35cm, push-pull valve. Non sterile. Pack of 10 bags. <i>(one pack of 10 pieces contains two straps with button)</i>	750ml	40 packs of 10 bags	
103404	Urine collection leg bag 750ml capacity with tube Ø 4.8x6.8mm, length 10cm, T-tap valve. Non sterile. Pack of 10 bags. <i>(one pack of 10 pieces contains two straps with button)</i>	750ml	40 packs of 10 bags	
103444	Urine collection leg bag 750ml capacity with tube Ø 4.8x6.8mm, length 35cm, T-tap valve. Non sterile. Pack of 10 bags. <i>(one pack of 10 pieces contains two straps with button)</i>	750ml	40 packs of 10 bags	
113411C1	Urine collection leg bag 750ml capacity with tube Ø 4.8x6.8mm, length 10cm, T-tap valve. Sterile.	750ml	100 bags individually packaged	
113448C1	Urine collection leg bag 750ml capacity with tube Ø 4.8x6.8mm, length 35cm, T-tap valve. Sterile.	750ml	100 bags individually packaged	
LC-SET2	Pack of n. 2 elastic straps to secure the bag to the leg, with holes to adjust the fitting. Equipped with closure button. Non sterile.	/	25 packs	

Euro 5 bag


5000ml T-tap valve

REF	Description	Capacity	Sales box	Picture
114541	Urine collection bag 5000ml capacity with connecting tube Ø 6x9mm, length 130cm, open jaw blue slide clamp, bag hanger, Pasteur drip chamber , hydrophobic filter, needle-free sampling port with protective cap, T-tap drain valve. Sterile.	5000ml	50 bags	
114543	Urine collection bag 5000ml capacity with connecting tube Ø 6x9mm, length 130cm, open jaw blue slide clamp, bag hanger, without Pasteur drip chamber , hydrophobic filter, needle-free sampling port with protective cap, T-tap drain valve. Sterile.	5000ml	50 bags	

Euro 5 emptying bag 5000ml without outlet

REF	Description	Capacity	Sales box	Picture
104561	Liquid emptying bag 5000ml capacity with connecting tube Ø 6x9mm, length 10cm, ON/OFF tube clamp, hydrophobic filter, without outlet, with soft female connector for connection to the main bag. Non Sterile.	5000ml	100 bags	
104568	Liquid emptying bag 5000ml capacity with connecting tube Ø 6x9mm, length 10cm, ON/OFF tube clamp, hydrophobic filter, without outlet, with soft female connector for connection to the main bag. Equipped with bag hanger. Non Sterile.	5000ml	100 bags	

Sterile collection bags



Description








- PVC sterile collection bag, available in various capacities: 500, 2000, 5000 and 8000 ml
- Non return valve
- Connection tube length 90 cm, supplied with proper connector
- With or without “T” cross outlet
- Disposable, sterile, double packaged

Advantages

- Wide choice of connectors for multiple applications
- High quality “T” valve assures leakage free bag
- Non return valve for safe use
- Sterile double packaging for aseptic use

Code	Description	Capacity	Q.ty / box
MS/00055x		2000 ml	100
MS/00056x		500 ml	100
MS/00057x	Sterile collection bag, with non return valve, cross outlet. Connecting tube length 90 cm	500 ml	100
MS/00058x	Connector as described in the table below	5000 ml	100
MS/00059x		8000 ml	100

Note: replace “X” in the code with the connector type, as described in the table below

“x” code specification	Description		“x” code specification	Description	
G	Big conic connector		P	Female Luer Lock	
M	Medium conic connector		U	Universal connector	
D	Small conic connector		LM	Male Luer Lock with cap	
L	Big Luer Lock				

CE 0373