

Sterile collection bags



Description








- PVC sterile collection bag, available in various capacities: 500, 2000, 5000 and 8000 ml
- Non return valve
- Connection tube length 90 cm, supplied with proper connector
- With or without "T" cross outlet
- Disposable, sterile, double packaged

Advantages

- Wide choice of connectors for multiple applications
- High quality "T" valve assures leakage free bag
- Non return valve for safe use
- Sterile double packaging for aseptic use

Code	Description	Capacity	Q.ty / box
MS/00055x		2000 ml	100
MS/00056x		500 ml	100
MS/00057x	Sterile collection bag, with non return valve, cross outlet. Connecting tube length 90 cm	500 ml	100
MS/00058x	Connector as described in the table below	5000 ml	100
MS/00059x		8000 ml	100

Note: replace "X" in the code with the connector type, as described in the table below

"x" code specification	Description		"x" code specification	Description	
G	Big conic connector		P	Female Luer Lock	
M	Medium conic connector		U	Universal connector	
D	Small conic connector		LM	Male Luer Lock with cap	
L	Big Luer Lock				

CE 0373