


Thoracentesis Paracentesis Basic Set



Description

- Thoracentesis - paracentesis set
- Provided with 50/60ml aspiration syringe and three-way stopcock
- Cannula needle length 50 or 80mm, size 14, 16 or 18G
- Transparent and graduated collection bag 2000ml capacity, with non-return valve and connection tube length 90cm
- Disposable double packaged sterile

Use

- Drainage of pleural (thoracentesis) or abdominal (paracentesis) fluids

Advantages

- Each set includes all the accessories needed for immediate use
- Wide range of needles available
- Collection bag with non-return valve for safe use
- Sterile double packaging for aseptic use

| Code | Description | Needle length | Bag capacity | Q.ty/box |
|--------------|----------------------------|---------------|--------------|----------|
| MS/310/50/14 | Set with 1 needle size G14 | 50mm | 2000ml | 80 |
| MS/310/50/16 | Set with 1 needle size G16 | 50mm | 2000ml | 80 |
| MS/310/50/18 | Set with 1 needle size G18 | 50mm | 2000ml | 80 |
| MS/310/80/14 | Set with 1 needle size G14 | 80mm | 2000ml | 80 |
| MS/310/80/16 | Set with 1 needle size G16 | 80mm | 2000ml | 80 |
| MS/310/80/18 | Set with 1 needle size G18 | 80mm | 2000ml | 80 |

