


Thoracentesis Paracentesis set



THORACENTESIS / PARACENTESIS SET WITH 3-WAY STOPCOCK



Three standard needles



Verres needle

- Each set includes all the accessories needed for immediate use
- Wide range of needles available
- Collection bag with cross outlet and non return valve
- Disposable, double packaged sterile



Cannula needle

Code	Description	Needle length	Bag capacity	Q.ty/box
MS/300/50	Set with 3 standard needles sizes G14, G16, G18	50 mm	2000 ml	80
MS/300/80	Set with 3 standard needles sizes G14, G16, G18	80 mm	2000 ml	80
MS/301/G	Set with 1 Verres needle size G14, G15 or G16	100 mm	2000 ml	50
MS/306/size	Set with 1 cannula needle (various sizes/lengths)	----	2000 ml	50
MS/302/50	Set with 3 standard needles sizes G14, G16, G18	50 mm	5000 ml	80
MS/302/80	Set with 3 standard needles sizes G14, G16, G18	80 mm	5000 ml	80
MS/303/G	Set with 1 Verres needle size G14, G15 or G16	100 mm	5000 ml	50
MS/307/size	Set with 1 cannula needle (various sizes/lengths)	----	5000 ml	50
MS/304/50	Set with 3 standard needles sizes G14, G16, G18	50 mm	8000 ml	80
MS/304/80	Set with 3 standard needles sizes G14, G16, G18	80 mm	8000 ml	80
MS/305/G	Set with 1 Verres needle size G14, G15 or G16	100 mm	8000 ml	50
MS/308/size	Set with 1 cannula needle (various sizes/lengths)	----	8000 ml	50

THORACENTESIS / PARACENTESIS SET WITH TWO NON RETURN VALVES

- Easy and fast operation thanks to two non return valves
- Ready to use set with wide choice of needles
- Collection bag with cross outlet and non return valve
- Disposable, double packaged, sterile




Code	Description	Needle length	Bag capacity	Q.ty/box
MS/300V/50	Set with 3 standard needles sizes G14, G16, G18	50 mm	2000 ml	80
MS/300V/80	Set with 3 standard needles sizes G14, G16, G18	80 mm	2000 ml	80
MS/301V/G	Set with 1 Verres needle size G14, G15 or G16	100 mm	2000 ml	40
MS/306V/size	Set with 1 cannula needle (various sizes/lengths)	----	2000 ml	80
MS/302V/50	Set with 3 standard needles sizes G14, G16, G18	50 mm	5000 ml	80
MS/302V/80	Set with 3 standard needles sizes G14, G16, G18	80 mm	5000 ml	80
MS/303V/G	Set with 1 Verres needle size G14, G15 or G16	100 mm	5000 ml	40
MS/307V/size	Set with 1 cannula needle (various sizes/lengths)	----	5000 ml	80
MS/304V/50	Set with 3 standard needles sizes G14, G16, G18	50 mm	8000 ml	80
MS/304V/80	Set with 3 standard needles sizes G14, G16, G18	80 mm	8000 ml	80
MS/305V/G	Set with 1 Verres needle size G14, G15 or G16	100 mm	8000 ml	40
MS/308V/size	Set with 1 cannula needle (various sizes/lengths)	----	8000 ml	40

SPARE COLLECTION BAG

- Wide choice of bag capacities
- Each bag includes a non return valve and a cross outlet
- 90 cm connecting tube
- Disposable, double packaged, sterile





2000 ml



5000 ml



8000 ml

Code	Description	Bag capacity	Connector	Q,ty/box
MS/00055P	Spare sterile collection bag, for set with 3-way stopcock, with female Luer Lock connector	2000 ml		100
MS/00058P		5000 ml		
MS/00059P		8000 ml		
MS/00055L	Spare sterile collection bag, for set with non return valves, with female Luer Lock connector	2000 ml		100
MS/00058L		5000 ml		
MS/00059L		8000 ml		

CE 0373