

Bracci ureteral catheter



Description

- PVC Bracci catheter for ureterocutaneostomy
- Length cm. 150, ORX
- With 9 holes at distal end
- Closed tip
- Available sizes: Ch. 06, 08, 10, 12,
- Disposable, sterile, double packaged

Use

- Urine drainage (ureterocutaneostomy)

Advantages

- Easy positioning thanks to the ORX line
- PVC hardness for use without internal guide
- Wide range of available sizes
- Sterile double packaging for aseptic use

Code	Description	Holes	Size	Length	Q.ty / box
MS/00075/06 R	Ureteral Bracci catheter	9	Ch. 06	150 cm	50
MS/00075/08 R			Ch. 08		
MS/00075/10 R			Ch. 10		
MS/00075/12 R			Ch. 12		

